


Efficacy Characteristics of Biocidal ZF™ (ready to use disinfectant spray)

Efficacy Characteristics of AquaClean in 1:200 aqueous dilution



[min₅]

minimal killing concentration
of active substance in suspension (H₂O)

5 minutes incubation time

0,03% Albumin load

water hardness <5° dH

[min₂₀]

minimal killing concentration
of active substance in suspension (H₂O)

20 minutes incubation time

0,03% Albumin load

water hardness <5° dH

F

multiple of minimal killing concentration
of active ingredients present
in ready to use Biocidal ZF™
(displayed separately
for each test organism)

range of individual
multiple killing concentrations
in Biocidal ZF™
range of individual
multiple killing concentrations
in AquaClean in 1:200 aqueous dilution

